



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## **Gilson PIPETMAN® L P2L microliter pipette, plastic ejector, volume 0.2-2 µL**

**€327.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Single-channel pipette, variable type

More productivity more convenience

- ergonomic design
- light in weight
- 3 tip ejector positions - ideal for left- and right-handed users

GLP Features

- Volume lock/unlock mechanism
- unique 2D identification code
- traceable serial number
- color coded push buttons

## Additional Information

No.	CF-0607
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4425
EAN	4058072222796
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	2 µl
Volume MIN	0.2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code with dispensing function ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Color code	orange
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L P2L
Division	0.002 µl
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Measurement deviation systematic in %	1.5 %
Measurement deviation systematic in µl	0.03 µl
Measurement deviation random in %	2 %
Measurement deviation random in µl	0.04 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

