



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L P10L microliter pipette, plastic ejector, volume 1-10 µL

**€327.00
plus VAT &
Shipping**

Product Images



Description

Single-channel pipette, variable type

More productivity more convenience

- ergonomic design
- light in weight
- 3 tip ejector positions - ideal for left- and right-handed users

GLP Features

- Volume lock/unlock mechanism
- unique 2D identification code
- traceable serial number
- color coded push buttons

Additional Information

No.	CF-0613
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4426
EAN	4058072222857
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	10 µl
Volume MIN	1 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code with dispensing function ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Color code	red
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L P10L
Division	0.01 µl
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.1 µl
Measurement deviation random in %	0.4 %
Measurement deviation random in µl	0.04 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

