



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L P10L microliter pipette, plastic ejector, volume 1-10 mL

€327.00
plus VAT &
Shipping

Product Images



Description

Single-channel pipette, variable type

More productivity more convenience

- ergonomic design
- light in weight
- 3 tip ejector positions - ideal for left- and right-handed users

GLP Features

- Volume lock/unlock mechanism
- unique 2D identification code
- traceable serial number
- color coded push buttons

Additional Information

No.	CF-0643
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4432
EAN	4058072223151
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	10000 µl
Volume MIN	1000 µl
Microliter pipettes features	with certification with GLP functions with color code with dispensing function ergonomically shaped IVD compliant
Number of channels	1
Operation via	Knobs/Buttons
Designation tip (Gilson)	D10ml
Energy supply	without power supply
Color code	light blue
Functionality	mechanical
Manufacturer	Gilson
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	60 µl
Measurement deviation random in %	0.16 %
Measurement deviation random in µl	16 µl
Model Microliter pipette/ Dispenser	Pipetman L P10mLL
Division	10 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

