



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Pipetman® electronic, P20M 2-20 µL

Special Price

€205.00 was

€270.00

**plus VAT &
Shipping**

Product Images



Description

Gilson Pipetman® M electronic

As easy to use as a mechanical pipette

- One-click operation like a mechanical pipette
- intuitive software for both experienced and inexperienced users
- user-independent pipetting

Pipetting modes: pipetting, aliquoting, mixing, reverse pipetting

GLP functions: Lockable volume, push button with color coding, service reminder

Additional Information

No.	CF-0882
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4434
EAN	4066292499263
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	2 µl
Volume MIN	20 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with mains/charging adapter with GLP functions with color code batteries/battery included IVD compliant
Number of channels	1
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Energy supply	Battery
Color code	yellow
Functionality	electric
Manufacturer	Gilson
Measurement deviation systematic in %	0.75 %
Measurement deviation systematic in µl	0.15 µl
Measurement deviation random in %	0.25 %
Measurement deviation random in µl	0.05 µl
Model Microliter pipette/ Dispenser	Pipetman M P20M
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

