



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Socorex Acura® electro 926 XS Electronic micropipette, volume 5-100 µl

**€604.24**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Socorex Acura® electro 926 XS electronic **microliter pipettes** are pipettes that make pipetting more flexible and safer than ever.

The electronic pipettes feature intuitive software that provides quick access to six work to six work functions and the different volume ranges.

Optimal ergonomics, low weight.

- Reduced length and tapered end of the shaft
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for the most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user

Two-year warranty.

## Additional Information

No.	SR-0009
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7205
EAN	4058072365615
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	100 µl
Volume MIN	5 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	1
Energy supply	Power grid
Functionality	electric
Manufacturer	Socorex
Scope of delivery	A pipette, individual control certificate, an instruction manual, samples of Qualitips® pipette tips
Measurement deviation systematic in %	0.5 %
Measurement deviation random in %	0.2 %
Model Microliter pipette/ Dispenser	Socorex Acura electro XS 926
Division	0.5 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	TÜV certificate

