



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® electro 936 Electronic macropipette, volume 0.25-5 ml

€628.68
plus VAT &
Shipping

Product Images



Description

The **Socorex Acura® electro 936 macropipettes** are electronic pipettes that make pipetting more flexible and safer than ever.

The electronic pipettes feature intuitive software that allows quick access to six work functions and the different volume ranges.

- Optimal ergonomics, low weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for the most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user
- Replaceable protection filter
- Adapter for Pasteur pipettes (2 ml) and straw tips

Two years warranty.

Additional Information

No.	SR-0021
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7217
EAN	4058072365721
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	5000 µl
Volume MIN	250 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with mains/charging adapter with charging station with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	1
Energy supply	Battery
Functionality	electric
Model Microliter pipette/ Dispenser	Socorex Acura electro XS 936
Division	0.05 ml
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	0.5 %
Measurement deviation systematic in µl	25 µl
Measurement deviation random in %	0.15 %
Measurement deviation random in µl	7.5 µl
Type pipette	Microliter pipettes
Type certificate	TÜV certificate

