



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® dilute 810 pipette for 1:10 dilutions, 1 / 0.1 ml

**€270.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Micropipette Acura® dilute for 1:10 dilutions**

The air cushion pipette with two pre-calibrated stages allows sequential dispensing of 1 and 0.1 ml of the same liquid. A universal metal cone ensures perfect fit of long PP straw tips, advantageous for sampling in narrow vessels (e.g. stomacher bags).

An ideal alternative to conventional graduated pipettes for standard 1:10 dilutions in bacteriology.

Three-year warranty.

- Two pre-calibrated fixed volumes - no adjustment necessary
- Smooth activation, excellent ergonomics
- Replaceable PE protection filter
- Justip™ system\* for tip ejection height adjustment
- Independent calibration of each volume
- Easy maintenance, cleaning and disinfection
- Fully autoclavable at 121°C

\* Socorex patented

## Additional Information

No.	SR-0114
Manufacturer (Brand)	Socorex
old neolab article no.	7-7363
EAN	4058072366445
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code with dispensing function ergonomically shaped
Number of channels	1
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	0.5 %
Measurement deviation systematic in µl	5 µl
Measurement deviation random in %	2.5 %
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

