



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® manual 835 Laboratory macropipette, volume 0.5-5 ml

**€331.76**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Socorex Acura® manual 835 macropipettes** were developed with the highest customer demands in mind. Suitable for indoor and outdoor environmental analysis applications. As well as for cell culture or clinical chemistry. They offer high ease of use and reliable results. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable to a very fine volume resolution. The extremely smooth volume adjustment allows simultaneous one-hand operation. The drop is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure the utmost stability of the volume setting within the pipetting process.

For the Acura® manual 835 models with 2 or 5 ml volume, a version including a Pasteur pipette adapter ( 7-7314/7-7316) is also available.

This is a useful addition to the use of PP tips for the use of Pasteur pipettes.

In addition, the pipette has a replaceable protective filter. This prevents liquid from entering the pipette.

- Excellent ergonomics
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection\* for the most common tips
- Integrated In-Lab calibration system Swift-Set\*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

\* Socorex patented

## Additional Information

No.	SR-0103
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7316
EAN	4058072366339
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	5000 µl
Volume MIN	500 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	1.5 %
Measurement deviation systematic in µl	75 µl
Measurement deviation random in %	1.5 %
Measurement deviation random in µl	75 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

