



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® manual 825 Laboratory micropipette, volume 0.5-10 µl

**€288.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Socorex Acura® manual microliter pipettes** have been developed with the highest customer demands in mind and guarantee convenience and reliable results for both routine and complicated pipetting processes. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable up to a very fine volume resolution. The extremely smooth volume adjustment allows simultaneous one-hand operation. The drop is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure extreme stability of volume adjustment within pipetting operations.

- Excellent ergonomics
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection\* for the most common tips
- Integrated In-Lab calibration system Swift-Set\*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

\* Socorex patented

## Additional Information

No.	SR-0080
Manufacturer (Brand)	Socorex
old neolab article no.	7-7282
EAN	4058072162832
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Manufacturer	Socorex
Measurement deviation systematic in %	2.5 %
Measurement deviation systematic in µl	0.25 µl
Measurement deviation random in %	1.8 %
Measurement deviation random in µl	0.18 µl
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

