



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® manual 825 Laboratory micropipette, volume 10-100 µl

€270.40
plus VAT &
Shipping

Product Images



Description

Socorex Acura® manual microliter pipettes have been developed with the highest customer demands in mind and guarantee convenience and reliable results for both routine and complicated pipetting processes. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable up to a very fine volume resolution. The extremely smooth volume adjustment allows simultaneous one-hand operation. The drop is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure extreme stability of volume adjustment within pipetting operations.

- Excellent ergonomics
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection* for the most common tips
- Integrated In-Lab calibration system Swift-Set*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

* Socorex patented

Additional Information

No.	SR-0084
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7286
EAN	4058072162870
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	100 µl
Volume MIN	10 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	1.5 %
Measurement deviation systematic in µl	1.5 µl
Measurement deviation random in %	1 %
Measurement deviation random in µl	1 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

