



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® manual 826 XS Laboratory micropipette, volume 0.5-10 µl

€319.80
plus VAT &
Shipping

Product Images



Description

Socorex Acura® manual microliter pipettes have been developed with the highest customer demands in mind and guarantee convenience and reliable results in both routine and complicated pipetting processes. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable up to a very fine volume resolution.

The **Acura manual 826XS models** are ideally **suited for research laboratories**. The short shaft allows for optimized instrument guidance. Ideal for precise applications, such as pipetting into microtubes. The innovative sealing ring concept* stands for extremely smooth pipetting and reduced hand fatigue during workflows. The extremely smooth volume adjustment allows simultaneous one-hand operation. The ejector is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure extreme stability of the volume setting during pipetting operations.

- Excellent ergonomics
- Reduced weight
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection* for the most common tips
- Integrated In-Lab calibration system Swift-Set*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

* Socorex patented

Additional Information

No.	SR-0089
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7300
EAN	4058072162924
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.1 µl
Measurement deviation random in %	0.5 %
Measurement deviation random in µl	0.05 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

