



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Socorex Acura® manual 826 XS Laboratory micropipette, volume 10-100 µl

**€300.04**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

**Socorex Acura® manual microliter pipettes** have been developed with the highest customer demands in mind and guarantee convenience and reliable results in both routine and complicated pipetting processes. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable up to a very fine volume resolution.

The **Acura manual 826XS models** are ideally **suited for research laboratories**. The short shaft allows for optimized instrument guidance. Ideal for precise applications, such as pipetting into microtubes. The innovative sealing ring concept\* stands for extremely smooth pipetting and reduced hand fatigue during workflows.

The extremely smooth volume adjustment allows simultaneous one-hand operation. The ejector is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure extreme stability of the volume setting during pipetting operations.

- Excellent ergonomics
- Reduced weight
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection\* for the most common tips
- Integrated In-Lab calibration system Swift-Set\*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

\* Socorex patented

## Additional Information

No.	SR-0093
Manufacturer (Brand)	Socorex
old neolab article no.	7-7304
EAN	4058072162962
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Volume MAX	100 µl
Volume MIN	10 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura manual
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.8 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	0.2 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

