



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® self-refill 865 Self-filling microdispensing pipette, volume 5-50 µl

€424.32
plus VAT &
Shipping

Product Images



Description

The **Acura® self-refill micro dispensing pipette** of the Socorex brand is ideally suited for serial dispensing in the microliter range. A three-way valve allows continuous dispensing. Carefully selected materials guarantee long service life and optimum resistance to aggressive media.

- User-friendly, extremely smooth operation of all functions
- Large volume display window
- Versatile liquid supply from bottle, tube or syringe
- Simple calibration system (Socorex patented)
- Individual serial number and quality control certificate
- Fully assembled autoclavable at 121° C
- 3 years warranty

Additional Information

No.	SR-0001
Manufacturer (Brand)	Socorex
old neoLab article no.	7-0255
EAN	4058072365554
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	light gray
Volume MAX	50 µl
Volume MIN	5 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with dispensing function ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Socorex Acura self-refill 5 - 50 µl
Division	0.1 µl
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	1.5 %
Measurement deviation systematic in µl	0.75 µl
Measurement deviation random in %	0.4 %
Measurement deviation random in µl	0.2 µl
Type pipette	Micro Dispenser Microliter pipettes

