



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® self-refill 865 Self-filling microdispensing pipette, volume 20-200 µl

**€388.96**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **Acura® self-refill micro dispensing pipette** of the Socorex brand is ideally suited for serial dispensing in the microliter range. A three-way valve allows continuous dispensing. Carefully selected materials guarantee long service life and optimum resistance to aggressive media.

- User-friendly, extremely smooth operation of all functions
- Large volume display window
- Versatile liquid supply from bottle, tube or syringe
- Simple calibration system (Socorex patented)
- Individual serial number and quality control certificate
- Fully assembled autoclavable at 121° C
- 3 years warranty

## Additional Information

No.	SR-0002
Manufacturer (Brand)	Socorex
old neolab article no.	7-0256
EAN	4058072365561
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	light gray
Volume MAX	200 µl
Volume MIN	20 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with dispensing function ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Manufacturer	Socorex
Measurement deviation systematic in %	1.8 %
Measurement deviation systematic in µl	3.6 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.6 µl
Model Microliter pipette/ Dispenser	Socorex Acura self-refill 20 - 200 µl
Division	0.2 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Micro Dispenser Microliter pipettes

