



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® electro 936 Electronic macropipette, starter pack, volume 0.1-2 ml

€904.80
plus VAT &
Shipping

Product Images



Description

The **Socorex Acura® electro 936 macropipettes** are flexible electronic pipettes that make pipetting more flexible and safer than ever.

The Acura® electro Initial Packs include an electronic pipette, an individual control certificate, operating instructions, a charging station, batteries with adapter and cable, and samples of Qualitips® pipette tips.

The electronic pipettes feature intuitive software that provides quick access to the six work to the six working functions and the different volume ranges.

Two-year warranty.

- Optimal ergonomics, light weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user
- Replaceable protection filter
- Adapter for Pasteur pipettes and stalk tips

Additional Information

No.	SR-0023
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7219
EAN	4058072365745
Transport temperature	Room temperature
autoclavable at 121°C	No
Volume MAX	2000 µl
Volume MIN	100 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with mains/charging adapter with charging station with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	1
Energy supply	Battery
Functionality	electric
Model Microliter pipette/ Dispenser	Socorex Acura electro XS 936
Division	0.01 ml
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Scope of delivery	An electronic pipette, an individual control certificate, operating instructions, a charging station, batteries with adapter and cable, and samples of Qualitips® pipette tips.
Measurement deviation systematic in %	0.5 %
Measurement deviation systematic in µl	10 µl
Measurement deviation random in %	0.15 %
Measurement deviation random in µl	3 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

