



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## sunlab® microliter pipette, mechanical, 2-20 $\mu$ l

**€135.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Sunlab® mechanical microliter pipette is the right pipette when it comes to fast and easy pipetting. It is lightweight, ergonomically designed and can be optimally used for all common tasks. The microliter pipette is fully autoclavable and available in different volumes from 0.1 to 10000 µl.

The volume of the pipette can be easily adjusted via a rotary wheel on the pipetting knob and read on the display. The ejector makes it easy to remove used pipette tips. The Sunlab® microliter pipette is calibrated according to the DIN EN ISO 8655 standard

## Additional Information

---

No.	D-8703
Manufacturer (Brand)	sunlab
EAN	4058072085834
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	20 µl
Volume MIN	2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Sunlab 1-channel microliter pipette
Division	0.5 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Sunlab
Sunlab® Model number	Microliter pipette SUI703
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

