



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific FinnpiPETTE™ F2 microliter pipette, volume 1-10 µl

**€343.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Single-channel pipette, variable type

The new Finnpiquette F2 combines all the features that are important for a pipette - easy operation, simple adjustment, complete autoclavability as well as exceptional reliability.

- novel volume setting mechanism for highest accuracy, precision and durability
- easy pipetting for fatigue-free work over a longer period of time
- multifunctional pipetting knob prevents unintentional volume adjustment during pipetting
- very large Ergovisio display - easy to read for easier volume adjustment
- completely autoclavable
- 5-year warranty when registered on the Internet

## Additional Information

No.	PO-0099
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1452
EAN	4058072338244
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	pink
Volume MAX	10 µl
Volume MIN	1 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Finnpipette F2 10µl micro
Division	0.01 µl
Designation tip (thermal)	Flex 10 50-micro
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Thermo Scientific
Measurement deviation systematic in %	1.2 %
Measurement deviation systematic in µl	0.12 µl
Measurement deviation random in %	0.8 %
Measurement deviation random in µl	0.08 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

