



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific FinnpiPETTE™ F2 microliter pipette, volume 2-20 µl

€343.00
plus VAT &
Shipping

Product Images



Description

Single-channel pipette, variable type

The new Finnpipette F2 combines all the features that are important for a pipette - easy operation, simple adjustment, complete autoclavability as well as exceptional reliability.

- novel volume setting mechanism for highest accuracy, precision and durability
- easy pipetting for fatigue-free work over a longer period of time
- multifunctional pipetting knob prevents unintentional volume adjustment during pipetting
- very large Ergovisio display - easy to read for easier volume adjustment
- completely autoclavable
- 5-year warranty when registered on the Internet

Additional Information

No.	PO-0101
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1453
EAN	4058072338251
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	turquoise
Volume MAX	20 µl
Volume MIN	2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Designation tip (thermal)	50-micro
Energy supply	without power supply
Functionality	mechanical
Manufacturer	Thermo Scientific
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	0.5 %
Measurement deviation random in µl	0.1 µl
Model Microliter pipette/ Dispenser	Finnpipette F1 20µl micro
Division	0.02 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

