



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific FinnpiPETTE F2 single channel, 2-20 $\mu$ l

**€343.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Single-channel pipette, variable type

The new Finnpipette F2 combines all the features that are important for a pipette - easy operation, simple adjustment, complete autoclavability as well as exceptional reliability.

- novel volume setting mechanism for highest accuracy, precision and durability
- easy pipetting for fatigue-free work over a longer period of time
- multifunctional pipetting knob prevents unintentional volume adjustment during pipetting
- very large Ergovisio display - easy to read for easier volume adjustment
- completely autoclavable
- 5-year warranty when registered on the Internet

## Additional Information

No.	PO-0102
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1454
EAN	4058072338268
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Volume MAX	20 µl
Volume MIN	2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Finnpipette F2 20µl
Division	0.02 µl
Designation tip (thermal)	Flex 200
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Thermo Scientific
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	0.5 %
Measurement deviation random in µl	0.1 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

