



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Transferpette® S, 20–200 µl, Variable, D-200, DE-M

Special Price

€185.57 was

€265.10

**plus VAT &
Shipping**

Product Images



Description

With the Transferpette® S piston-stroke pipette, you can work efficiently and ergonomically with small and large volumes. Low operating forces and a short stroke path ensure relaxed pipetting even with long pipetting series. Thanks to its perfect ergonomics, the Transferpette® S fits comfortably in the hand in any holding position, whether with the right or left hand, whether with large or small hands. There is no need to reach around, because the microliter pipette offers true one-handed operation:

The volume adjustment protection is easy to operate even when wearing gloves and secures your set volume against accidental adjustment. To make adjustment easy, the Transferpette® S also features Easy Calibration technology.

Adjustment can thus be carried out without tools. And to ensure that you always order the right pipette tips, the compatibility table shows you which tips fit which instrument. At a glance-

Volume ranges single-channel pipette:

from 0.1 µl to 10 ml- Volume ranges multichannel pipette:
from 0.5 µl to 300 µl- Short stroke and low operating forces- Universal tip mounting cone- Fully autoclavable at 121 °C- Single-handed operation for right- and left-handed users- Color code on volume adjustment wheel, Shaft and volume adjustment protection- Easy tip attachment and ejection thanks to O-rings on multi-channel pipettes- Easy to adjust without tools thanks to Easy Calibration technology- Slimmer shaft for removal from narrow vessels- Operating temperature from +15 °C to +40 °C- Vapor pressure up to 500 mbar- Viscosity up to 260 mPa s

Additional Information

No.	BN-2659
Article - USP	In any posture: comfortable one-hand operation For all hands: ergonomic design Safe: volume adjustment protection prevents accidental adjustment
Manufacturer (Brand)	Brand
VGKL number	104975878
EAN	4058072217761
Gross weight	0.3025
HSNumber	84798997
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.18740 kg
Volume MAX	200 µl
Volume MIN	20 µl
Country of origin	Germany
Packing width	0.12000 m
Packing height	0.04800 m
Packing depth	0.28700 m
Packaging volume	1653.120000 ccm
Customs tariff number	84798997
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	1.2 µl
Measurement deviation random in µl	0.4 µl
Youtube URL	https://youtu.be/YLhNM-Cj5JQ

