



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoPette® microliter pipette 500–5000 µl

Special Price

€95.00 was

€142.50

**plus VAT &
Shipping**

Product Images



Description

Single-channel pipette, variable type

The new **neoPette** is a precision instrument of high quality with a distinctive design, covering a volume range from 0.5 to 10,000 μl . It features a continuously variable, digital volume adjustment.

Each pipette model is marked with its maximum volume on the handle (clip) for easy identification.

The neoPette bears the H conformity mark according to DIN EN ISO 8655-6.

Recalibration directly in the laboratory is easy thanks to Easy Calibration.

All models have a tip ejector. Completely autoclavable at 121°C.

Additional Information

No.	7-2135
Manufacturer (Brand)	neoPette
EAN	4066292552951
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Volume MAX	5000 µl
Volume MIN	500 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	neoPette digital
Division	5 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	neoLab
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	30 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	10 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 8655-2:2009-07

