



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoPette® microliter pipette 1-10 ml

Special Price

€95.00 was

€142.50

**plus VAT &
Shipping**

Product Images



Description

Single-channel pipette, variable type

The new **neoPette** is a precision instrument of high quality with a distinctive design, covering a volume range from 0.5 to 10,000 μl . It features a continuously variable, digital volume adjustment.

Each pipette model is marked with its maximum volume on the handle (clip) for easy identification.

The neoPette bears the H conformity mark according to DIN EN ISO 8655-6.

Recalibration directly in the laboratory is easy thanks to Easy Calibration.

All models have a tip ejector. Completely autoclavable at 121°C.

Additional Information

No.	7-2136
Manufacturer (Brand)	neoPette
EAN	4066292552968
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Volume MAX	10000 µl
Volume MIN	1000 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Manufacturer	neoLab
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	60 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	20 µl
Model Microliter pipette/ Dispenser	neoPette digital
Division	10 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 8655-2:2009-07

