



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## neoPette® Microliter pipette, electronic, 30–300 $\mu$ l

**€159.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

### **The revolution in liquid handling.**

As an electronic microliter pipette, the new neoPette® e offers everything you need for convenient pipetting. You can pipette, dispense, dilute and mix with it. Unlike conventional electronic microliter pipettes, the neoPette® e is shaped like a manual pipette and fits comfortably in the hand. Operation is simple and intuitive. Slim ergonomic design, low weight despite rechargeable battery. The neoPette® e is calibrated and adjusted according to the DIN EN ISO 8655 standard.

In addition to the neoPette® e, the scope of delivery includes a charging adapter and a charging stand for attachment to the shelf.

## Additional Information

No.	7-2142
Manufacturer (Brand)	neoPette
EAN	4058072083250
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Volume MAX	300 µl
Volume MIN	30 µl
sterile	No
Microliter pipettes features	with certification with mains/charging adapter with charging station including batteries/battery ergonomically shaped
Number of channels	1
Operation via	Knobs/Buttons
Energy supply	Battery
Functionality	electronic
Manufacturer	neoLab
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	1.8 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	0.6 µl
Model Microliter pipette/ Dispenser	neoPette digital
Division	2 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

