



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reference® 2, 12-channel, variable, 0.5 – 10 µL, medium gray

€1,135.00
plus VAT &
Shipping

Product Images



Description

The Reference® 2 microliter multipipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

- Simple and visual volume adjustment with volume fixation
- Volume range from 0.5 µl to 300 µl
- Fully autoclavable
- Single button operation
- reproducible seating of the tip on the pipette cone

Additional Information

No.	VB-1532
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0473
EAN	4058072416409
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	medium gray
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped
Number of channels	12
Operation via	Knobs/Buttons
Energy supply	without power supply
Color code	medium gray
Functionality	mechanical
Manufacturer	Eppendorf
Measurement deviation systematic in %	2 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	1 %
Measurement deviation random in µl	0.1 µl
Model Microliter pipette/ Dispenser	Reference 2 0.5 - 10µl
Division	0.01 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer

