



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reference® 2, 8-channel, variable, 10 - 100 µL, yellow

€1,037.00
plus VAT &
Shipping

Product Images



Description

The Reference® 2 microliter multipipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

- Simple and visual volume adjustment with volume fixation
- Volume range from 0.5 µl to 300 µl
- Fully autoclavable
- Single button operation
- reproducible seating of the tip on the pipette cone

Additional Information

No.	VB-1533
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0471
EAN	4058072416416
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Volume MAX	100 µl
Volume MIN	10 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped
Number of channels	8
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Model Microliter pipette/ Dispenser	Reference 2 10 - 100µl
Division	0.1 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.8 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.3 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer

