



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reference® 2, 12-channel, variable, 10 - 100 µL, yellow

€1,135.00
plus VAT &
Shipping

Product Images



Description

The Reference® 2 microliter multipipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

- Simple and visual volume adjustment with volume fixation
- Volume range from 0.5 µl to 300 µl
- Fully autoclavable
- Single button operation
- reproducible seating of the tip on the pipette cone

Additional Information

No.	VB-1534
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0474
EAN	4058072416423
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Volume MAX	100 µl
Volume MIN	10 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped
Number of channels	12
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Model Microliter pipette/ Dispenser	Reference 2 10 - 100µl
Division	0.1 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.8 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.3 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer

