



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Research® plus, 8-channel, variable; 0.5-10 µl, medium gray

Special Price  
**€705.60** was  
**€882.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Multichannel pipettes, type Variable

offer the same features as the <sup>Research®</sup> plus single-channel pipettes plus:

- Channel display: use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility for your individual needs

## Additional Information

No.	VB-0962
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0940
EAN	4058072170219
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	medium gray
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	8
Energy supply	without power supply
Color code	medium gray
Functionality	mechanical
Model Microliter pipette/ Dispenser	Research plus 8x10
Division	0.01 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	2 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	1 %
Measurement deviation random in µl	0.1 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

