



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Research® plus, 12-channel, variable; 0.5–10 µl, medium gray

Special Price

€772.00 was

€965.00

**plus VAT &
Shipping**

Product Images



Description

Multichannel pipettes, type Variable

offer the same features as the ^{Research®} plus single-channel pipettes plus:

- Channel display: use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility for your individual needs

Additional Information

No.	VB-0963
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0943
EAN	4058072169985
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	medium gray
Volume MAX	10 µl
Volume MIN	0.5 µL
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	12
Energy supply	without power supply
Color code	medium gray
Functionality	mechanical
Model Microliter pipette/ Dispenser	Research plus 12x10
Division	0.01 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	2 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	1 %
Measurement deviation random in µl	0.3 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

