



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Research® plus, 8-channel, variable; 30-300µl, orange

Special Price
€705.60 was
€882.00
plus VAT &
Shipping

Product Images



Description

Multichannel pipettes, type Variable

offer the same features as the ^{Research®} plus single-channel pipettes plus:

- Channel display: use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility for your individual needs

Additional Information

No.	VB-0966
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-0942
EAN	4058072411114
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	orange
Volume MAX	300 µl
Volume MIN	30 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	8
Energy supply	without power supply
Color code	orange
Functionality	mechanical
Model Microliter pipette/ Dispenser	Research plus 8x300
Division	0.2 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	1.8 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.9 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

