



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Research® plus, variable, 12- channel; 30-300 µl, orange

Special Price

**€772.00** was

**€965.00**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Multichannel pipettes, type Variable

offer the same features as the <sup>Research®</sup> plus single-channel pipettes plus:

- Channel display: use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility for your individual needs

## Additional Information

No.	VB-0967
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0945
EAN	4058072411121
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	orange
Volume MAX	300 µl
Volume MIN	30 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	12
Energy supply	without power supply
Color code	orange
Functionality	mechanical
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	1.8 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.9 µl
Model Microliter pipette/ Dispenser	Research plus 12x300
Division	0.2 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

