



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Research® plus, 8-channel, variable; 120-1200 µl, green

Special Price
€705.60 was
€882.00
plus VAT &
Shipping

Product Images



Description

Multichannel pipettes, type Variable

offer the same features as the ^{Research®} plus single-channel pipettes plus:

- Channel display: use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility for your individual needs

Additional Information

No.	VB-2497
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-0949
EAN	4058072425135
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1200 µl
Volume MIN	120 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	8
Energy supply	without power supply
Color code	green
Functionality	mechanical
Division	0.1 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

