



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson PIPETMAN® M P8X300M BT, 8-channel connected, Bluetooth®-enabled electronic pipette, volume 10–300 µL

**€1,129.00**  
**plus VAT &**  
**Shipping**

### Product Images



# Description

---

## **Electronic pipette PIPETMAN® M Connected, Bluetooth®-enabled.**

PIPETMAN® M Connected is a Bluetooth®-enabled programmable electronic pipette, as easy to use as a mechanical pipette.

- The electronic pipette with guaranteed performance, even in repetitive mode, according to ISO8655 specifications
- Programming of individual protocols (also possible via PC and PipettPilot™ app) for reliable pipetting. The motorized piston ensures high reproducibility.
- The intuitive user interface has five pipetting modes for numerous applications: "pipette", "repetitive", "mix", "reverse" and "custom".
- Aspirate and dispense with the click of a button and minimal effort. Combined with the ergonomic design, this provides maximum pipetting comfort, reducing pipetting fatigue and increasing pipetting efficiency.

GLP features:

- Lockable volume
- Push button with color coding
- Service reminder

## Additional Information

No.	CF-0806
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4284
EAN	4058072224714
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	300 µl
Volume MIN	10 µl
Operation via	Knobs/Buttons
Number of channels	8
Energy supply	Battery
Color code	green
Functionality	electric
Designation tip (Gilson)	D200 D300
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Scope of delivery	<ul style="list-style-type: none"><li>- PIPETMAN M Connected pipette</li><li>- Power Supply 5V with AC adapter and cable</li><li>- 4 Battery side tags</li><li>- 4 Ejector side tags</li><li>- Tip ejector extension (P10M model only)</li><li>- Quick Guide</li></ul>
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

