



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson PIPETMAN® G P12x10G, 12-channel microliter pipette, metal ejector, volume 1-10 µL

**€915.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

**multichannel Pipetman® G** Premium comfort in Gilson quality

- Minimal force required for pipetting
- Reduction of tip ejection forces
- Color coded push button
- Reliable pipetting from as little as 10% of nominal volume

## Additional Information

No.	CF-0782
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4261
EAN	4058072224486
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	10 µl
Volume MIN	1 µl
Operation via	Knobs/Buttons
Number of channels	12
Energy supply	without power supply
Color code	red
Functionality	mechanical
Designation tip (Gilson)	D10 DL10
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Scope of delivery	<p>pipETMAN G</p> <ul style="list-style-type: none"><li>● PIPETMAN G Quick Guide</li><li>● Safety bag</li><li>● Ejector spacers for D10 tips (for Px10G models only, refer to Ejector Spacers for PIPETMAN G Multichannel 10µl)</li></ul> certificate of Conformity (including barcode sticker)
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

