



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G P12x300G, 12-channel microliter pipette, metal ejector, volume 30–300 µL

€915.00
plus VAT &
Shipping

Product Images



Description

multichannel Pipetman® G Premium comfort in Gilson quality

- Minimal force required for pipetting
- Reduction of tip ejection forces
- Color coded push button
- Reliable pipetting from as little as 10% of nominal volume

Additional Information

No.	CF-0788
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4267
EAN	4058072224530
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	300 µl
Volume MIN	30 µl
Operation via	Knobs/Buttons
Number of channels	12
Energy supply	without power supply
Color code	green
Functionality	mechanical
Designation tip (Gilson)	D200 D300
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Scope of delivery	<p>pipETMAN G</p> <ul style="list-style-type: none">● PIPETMAN G Quick Guide● Safety bag● Ejector spacers for D10 tips (for Px10G models only, refer to Ejector Spacers for PIPETMAN G Multichannel 10µl) <p>certificate of Conformity (including barcode sticker)</p>
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

