



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G P8x10G, 8-channel microliter pipette, metal ejector, volume 1-10 µL

€777.00
plus VAT &
Shipping

Product Images



Description

multichannel Pipetman® G Premium comfort in Gilson quality

- Minimal force required for pipetting
- Reduction of tip ejection forces
- Color coded push button
- Reliable pipetting from as little as 10% of nominal volume

Additional Information

No.	CF-0781
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4260
EAN	4058072224479
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	10 µl
Volume MIN	1 µl
sterile	No
Length	275 mm
Operation via	Knobs/Buttons
Number of channels	8
Energy supply	without power supply
Color code	red
Functionality	mechanical
Designation tip (Gilson)	D200 D300
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Scope of delivery	<p>pipETMAN G</p> <ul style="list-style-type: none">● PIPETMAN G Quick Guide● Safety bag● Ejector spacers for D10 tips (for Px10G models only, refer to Ejector Spacers for PIPETMAN G Multichannel 10µl) certificate of Conformity (including barcode sticker)
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

