



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G P8x300G, 8-channel microliter pipette, metal ejector, volume 30–300 µL

€777.00
plus VAT & Shipping

Product Images



Description

multichannel Pipetman® G Premium comfort in Gilson quality

- Minimal force required for pipetting
- Reduction of tip ejection forces
- Color coded push button
- Reliable pipetting from as little as 10% of nominal volume

Additional Information

| | |
|----------------------------------|---|
| No. | CF-0787 |
| Manufacturer (Brand) | Gilson |
| old neoLab article no. | 7-4266 |
| EAN | 4058072224523 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | blue |
| Material | Polyvinylidene fluoride (PVDF) |
| Volume MAX | 300 µl |
| Volume MIN | 30 µl |
| sterile | No |
| Length | 260 mm |
| Operation via | Knobs/Buttons |
| Number of channels | 8 |
| Energy supply | without power supply |
| Color code | green |
| Functionality | mechanical |
| Designation tip (Gilson) | D200 D300 |
| Type Microliter pipette (design) | Air cushion pipette |
| Manufacturer | Gilson |
| Scope of delivery | <p>pipETMAN G</p> <ul style="list-style-type: none">● PIPETMAN G Quick Guide● Safety bag● Ejector spacers for D10 tips (for Px10G models only, refer to Ejector Spacers for PIPETMAN G Multichannel 10µl) certificate of Conformity (including barcode sticker) |
| Type pipette | Microliter pipettes |
| Type certificate | Declaration of conformity (standardization) |

