



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L P8X20L, 8-channel microliter pipette, plastic ejector, volume 2-20 µL

**€905.00
plus VAT &
Shipping**

Product Images



Description

a patented volume locking system guarantees accuracy and precision during long, intensive pipetting with maximum comfort volume adjustment during pipetting can be excluded has an innovative tip ejection system that offers the lowest ejection forces using a lever arm for all tip models the tip ejector ejector button is adjustable so that both right-handed and left-handed users can comfortably eject the tips particularly light in weight

Additional Information

No.	CF-0646
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4452
EAN	4058072223182
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	20 µl
Volume MIN	2 µl
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	8
Operation via	Knobs/Buttons
Designation tip (Gilson)	D200 DF200ST
Energy supply	without power supply
Color code	white
Functionality	mechanical
Manufacturer	Gilson
Measurement deviation systematic in %	2 %
Measurement deviation systematic in µl	0.4 µl
Measurement deviation random in %	1 %
Measurement deviation random in µl	0.2 µl
Model Microliter pipette/ Dispenser	Pipetman L 8x20 L
Division	0.02 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

