



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Gilson PIPETMAN® L P8X10L, 8-channel
microliter pipette, plastic ejector, volume 0.5-10
µL**

**€905.00
plus VAT &
Shipping**

Product Images



Description

a patented volume locking system guarantees accuracy and precision during long, intensive pipetting with maximum comfort volume adjustment during pipetting can be excluded has an innovative tip ejection system that offers the lowest ejection forces using a lever arm for all tip models the tip ejector ejector button is adjustable so that both right-handed and left-handed users can comfortably eject the tips particularly light in weight

Additional Information

No.	CF-0658
Manufacturer (Brand)	Gilson
old neolab article no.	7-4450
EAN	4058072223304
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyvinylidene fluoride (PVDF)
Volume MAX	10 µl
Volume MIN	0.5 µL
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	8
Energy supply	without power supply
Color code	red
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L 8x10 L
Division	0.01 µl
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Measurement deviation systematic in %	2 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	2 %
Measurement deviation random in µl	0.2 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

