



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® electro 956 Electronic 12- channel microliter pipette, volume 2.5-50 µl

€1,088.88
plus VAT &
Shipping

Product Images



Description

The **Socorex Acura® electro 956 Multi multichannel pip**ettes are flexible electronic pipettes that make pipetting more flexible and safer than ever.

The electronic pipettes feature intuitive software that allows quick access to the six working functions and the different to the six working functions and the different volume ranges.

Two-year warranty.

- Optimal ergonomics, low weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user

Additional Information

No.	SR-0031
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7227
EAN	4058072365820
Transport temperature	Room temperature
autoclavable at 121°C	No
Volume MAX	50 µl
Volume MIN	2.5 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	12
Energy supply	Battery
Functionality	electric
Manufacturer	Socorex
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.4 µl
Measurement deviation random in %	0.4 %
Measurement deviation random in µl	0.2 µl
Model Microliter pipette/ Dispenser	Socorex Acura electro 956
Division	0.25 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

