



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® electro 956 Electronic 8- channel microliter pipette, volume 2.5-50 µl

€996.32
plus VAT &
Shipping

Product Images



Description

The **Socorex Acura® electro 956 Multi multichannel pip**ettes are flexible electronic pipettes that make pipetting safer than ever.

The electronic pipettes feature intuitive software that allows quick access to six work functions and the different volume ranges.

- Optimal ergonomics, low weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejection for most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user
- Multichannel pipette housings can be rotated through 360
- Modular system, convertible with single-channel volume units

Two-year warranty.

Additional Information

No.	SR-0027
Manufacturer (Brand)	Socorex
old neolab article no.	7-7223
EAN	4058072365783
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	50 µl
Volume MIN	2.5 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	8
Energy supply	Battery
Functionality	electric
Model Microliter pipette/ Dispenser	Socorex Acura electro 956
Division	0.25 µl
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Socorex
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.12 µl
Measurement deviation random in %	0.4 %
Measurement deviation random in µl	0.2 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

