



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Acura® electro 956 Electronic 12- channel microliter pipette, volume 20-350 µl

**€1,364.48**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Socorex Acura® electro multichannel pipettes 956 Multi** are flexible electronic pipettes that make pipetting more flexible and safer than ever.

The pipettes are equipped with intuitive software that allows quick access to the six to the six working functions and the different volume ranges

**Acura® electro Initial Packs** include an electronic pipette, an individual control certificate, an operating manual, a charging station, batteries with adapter and cable, and samples of Qualitips® pipette tips.

Two-year warranty.

- Optimal ergonomics, low weight
- Easy to understand, logical program sequences
- Large display, reversible right / left hand
- Quickly replaceable handle battery
- Fast charging, long working autonomy
- Adjustable ejector\* for most common tips
- Pipetting cycle counter
- Interchangeable volume units, easy conversion by the user

## Additional Information

No.	SR-0041
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7237
EAN	4058072365929
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	350 µl
Volume MIN	20 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with titration function with mains/charging adapter with charging station with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	12
Energy supply	Battery
Functionality	electric
Manufacturer	Socorex
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	2.1 µl
Measurement deviation random in %	0.25 %
Measurement deviation random in µl	0.875 µl
Model Microliter pipette/ Dispenser	Socorex Acura electro 956
Division	5 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

