



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**sunlab® microliter pipette, 8-channel,
mechanical, 5-50 µl**

**€387.30
plus VAT &
Shipping**

Product Images



Description

The Sunlab® 8-channel mechanical microliter pipette is the right pipette for fast and easy pipetting. It is lightweight, ergonomically designed and can be optimally used for all common tasks. The microliter pipette is fully autoclavable and available in different volumes.

The volume of the pipette can be easily adjusted via a rotary wheel on the pipetting knob and read on the display. The ejector makes it easy to remove used pipette tips. The Sunlab® microliter pipette is calibrated according to the DIN EN ISO 8655 standard.

Additional Information

No.	D-8711
Manufacturer (Brand)	sunlab
EAN	4058072086459
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	50 µl
Volume MIN	5 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification ergonomically shaped
Number of channels	8
Energy supply	without power supply
Functionality	mechanical
Manufacturer	Sunlab
Model Microliter pipette/ Dispenser	Sunlab 8-channel microliter pipette
Sunlab® Model number	8-channel microliter pipette SUI711
Division	0.5 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

