



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**sunlab® microliter pipette, 8-channel,
mechanical, 50-300 µl**

**€387.30
plus VAT &
Shipping**

Product Images



Description

The Sunlab® 8-channel mechanical microliter pipette is the right pipette for fast and easy pipetting. It is lightweight, ergonomically designed and can be optimally used for all common tasks. The microliter pipette is fully autoclavable and available in different volumes .

The volume of the pipette can be easily adjusted via a rotary wheel on the pipetting knob and read on the display. The ejector makes it easy to remove used pipette tips. The Sunlab® microliter pipette is calibrated according to the DIN EN ISO 8655 standard.

Additional Information

No.	D-8712
Manufacturer (Brand)	sunlab
EAN	4058072086466
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	300 µl
Volume MIN	50 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification ergonomically shaped
Number of channels	8
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	Sunlab 8-channel microliter pipette
Division	5 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Sunlab
Sunlab® Model number	8-channel microliter pipette SU1712
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from the manufacturer

