



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finnpiquette™ F1 multichannel pipette, volume 5-50 µl

€951.00
plus VAT &
Shipping

Product Images



Description

Finnpipette F1 multichannel pipette is available in 8-channel version for volume range 5–50 µl.

As with the Finnpipette F1 single-channel pipettes, a novel volume setting mechanism ensures a high degree of accuracy and precision.

The soft-touch tip ejector makes work easier by significantly reducing the force required to eject multiple tips.

The super-blow-out function ensures accurate dispensing even at the lowest volume for the small-volume models.

Additional Information

No.	PO-0176
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1441
EAN	4066292155701
Transport temperature	Room temperature
autoclavable at 121°C	No
Volume MAX	50 µl
Volume MIN	5 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification ergonomically shaped
Number of channels	8
Energy supply	without power supply
Functionality	mechanical
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Thermo Scientific
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

