



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific FinnpiPETTE™ F1 multichannel pipette, volume 1-10 µl

€1,077.00
plus VAT &
Shipping

Product Images



Description

Finnpipette FI multichannel pipette is available in 12-channel version. As with Finnpipette FI single-channel pipettes, a novel volume adjustment mechanism ensures a high degree of accuracy and precision .

The soft-touch tip ejector makes work easier by significantly reducing the force required to eject multiple tips.

The super-blow-out function ensures accurate dispensing even at the lowest volume for the small-volume models.

Additional Information

No.	PO-0179
Manufacturer (Brand)	Thermo Scientific
old neolab article no.	3-1443
EAN	4058072338329
Transport temperature	Room temperature
autoclavable at 121°C	No
Volume MAX	10 µl
Volume MIN	1 µl
Microliter pipettes features	with certification ergonomically shaped
Number of channels	12
Operation via	Knobs/Buttons
Energy supply	without power supply
Functionality	mechanical
Manufacturer	Thermo Scientific
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

