



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reference® 2, 3-pack, Option 2: 2 - 20 µL yellow, 20 - 200 µL, 100 - 1,000 µL

€882.00
plus VAT &
Shipping

Product Images



Description

The Reference® 2 microliter pipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

3 variable adjustable Eppendorf Reference® 2 per pack and one additional accessory selectable.

Eppendorf Reference® 2 ; 3-pack option 2

3 microliter pipettes (2 - 20 µl, 20 - 200 µl, 100 - 1000 µl) + 1 piece epT.I.P.S.® box with matching pipette tips + ballpoint pen

Additional Information

No.	VB-1512
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0477
EAN	4058072416201
Transport temperature	Room temperature
Microliter pipettes features	with certification IVD Conform
Functionality	mechanical
Model Microliter pipette/ Dispenser	Reference 2 2 - 20µl Reference 2 20 - 200µl Reference 2 100 - 1000µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization) IVD Directive 98/79/EC

