



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Multipette M4 Sarter Kit (incl. Combitip Rack, Combitip Assortment Pack 1 μ l-10 ml)

€437.00
plus VAT &
Shipping

Product Images



Description

With variable volume selection, you can choose from up to 100 volumes. Combitips advanced also allow dispensing in up to 100 increments with just one filling, depending on the volume. The Multipette®M4 covers volume sizes from 0.1 ml to 50 ml with 9 Combitips advanced.

A volume spectrum from 1 µl to 10 ml.

- simplified volume setting with volume selection wheel
- Volume display with Combitip recognition via display
- Direct displacement principle
- easy Combitip pick-up
- ergonomically shaped operating and filling lever
- Motion sensor switches off the Multipette during breaks in use

Dispenser Multipette®M4 Starter Kit consists of:

- 1 Combitip Rack
- 1 Combitip assortment pack 1 µl-10 ml

Additional Information

No.	VB-1691
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-4288
EAN	4058072417802
Transport temperature	Room temperature
Functionality	mechanical
Model Microliter pipette/ Dispenser	Multipette M4
Type Microliter pipette (design)	Direct displacement pipette

