



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G Starter Kit, with P2G, P20G, P200G, P1000G microliter pipettes (metal ejector) and pipette tips

€1,049.00
plus VAT &
Shipping

Product Images



Description

Cost-effective Starter Kit with 4 pipettes.

P2G, P20G, P200G, P1000G (metal ejector) and Diamond Tips DL10, D200, D1000.

Included in the kit are 4 Single Pipette Holders, 4 Comfort Handle+12 Stickers, 1 Gilson Guide to Pipetting, 1 Two-Minute Inspection Poster and 1 User Guide.

Additional Information

No.	CF-0407
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4209
EAN	4058072162696
Transport temperature	Room temperature
autoclavable at 121°C	No
Number of channels	1
Designation tip (Gilson)	DL10 D200 D1000
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman G P2G Pipetman G P20G Pipetman G P200G Pipetman G P1000G
Type Microliter pipette (design)	Air cushion pipette

