



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L Starter Kit, with P2L, P20L, P200L, P1000L Microliter pipettes (metal ejector) and pipette tips

€1,178.00
plus VAT &
Shipping

Product Images



Description

Low cost starter kit with 4 pipettes.

P2L, P20L, P200L, P1000L (metal ejector) and Diamond Tips DL10, D200, D1000.

Included in the kit are 4 single pipette holders, 1 Gilson Guide to Pipetting, 1 Two-Minute Inspection Poster and 1 User Guide.

Additional Information

No.	CF-0408
Manufacturer (Brand)	Gilson
old neolab article no.	7-4208
EAN	4058072162702
Transport temperature	Room temperature
autoclavable at 121°C	No
Number of channels	1
Designation tip (Gilson)	DL10 D200 D1000
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L P2L Pipetman L P20L Pipetman L P200L Pipetman L P1000L
Type Microliter pipette (design)	Air cushion pipette

