



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson PIPETMAN® G Starter Kit, with P2G, P10G, P100G microliter pipettes (metal ejector) and pipette tips

**€794.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Inexpensive in one set

### **Order No. 7-4424**

Pipetman G Microvolume Kit

- Gilson Pipetman G P2G, P10G and P100G
- Tips DL10 and D200 in Tippack Rack
- 3 single holders
- 3 plastic tip ejectors
- 1 Gilson Guide to Pipetting

## Additional Information

---

No.	CF-0414
Manufacturer (Brand)	Gilson
old neolab article no.	7-4424
EAN	4058072221867
Transport temperature	Room temperature
autoclavable at 121°C	No
Number of channels	1
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman G P2G Pipetman G P10G Pipetman G P100G
Type Microliter pipette (design)	Air cushion pipette

